

3 September 2022

File Memo

Subject: Maintaining the Zoll AED Plus Battery

Website: <https://www.aedusa.com/knowledge/zoll-aed-plus-battery-replacement-instructions/>



The battery should be replaced every 2-5 years. Uses Duracell Lithium 123 batteries (10). Check the date on the batteries to plan the next replacement.

<https://www.aedusa.com/product/zoll/zoll-aed-plus-batteries>



How to change the 10 batteries: VIDEO https://www.youtube.com/watch?v=xa5Ez_jlIfM&t=30s

Check also the expiration date on the pads. Order new pads when necessary. Be sure to use pads compatible with the PLUS model.

Steve Moese, President, Cincinnati Bridge Association.